

# Report

## Resolution of the planet



## Introduction

The **planet** is the high resolution portable X-ray diffractometer from xplorex. With ‘high resolution’ we mean that the attainable resolution with the planet is comparable to a standard laboratory system and thus by far the best resolution hitherto obtained with a portable system.

In this note we demonstrate the attainable resolution with the planet with a normal measurement on normal commercial LaB<sub>6</sub>.

## Experimental

About 1g of the as-received powder was pressed into the sample holder (so-called front loading). The prepared specimen was mounted in the diffractometer and measured. Because the **planet** uses a reflection geometry, sample preparation is not critical even for heavy absorbers like LaB<sub>6</sub>. The measurement parameters are summarized in Table 1.

Parameter	Value
<b>Range</b>	
Start Angle ( $^{\circ}2\theta$ )	19.3
End Angle ( $^{\circ}2\theta$ )	101.3
Step size	Variable; The <b>planet</b> has predefined step sizes
Angle of incidence ( $^{\circ}$ )	14.47
<b>Integration time</b>	120s / data point (the measurement comprises 14 times 640 data points)
<b>Total Measurement time</b>	32 minutes
<b>Diffractometer settings</b>	
High Tension	30 kV
Emission Current	0.65 mA
Tube anode	Cu
Focus dimensions	40 $\mu$ m diameter
Take off angle ( $^{\circ}$ )	8
Beam divergence ( $^{\circ}$ )	1.25
<b>Specimen dimensions</b>	7mm diameter; 2mm thickness
<b>Spinning frequency</b>	0.5 Hz
<b>Optical path</b>	Seemann – Bohlin based
<b>Focusing circle radius</b>	160mm
<b>Detector</b>	Dectris’ Mythen 1D solid state linear detector
<b>Identification software</b>	Match! From Crystal Impact
<b>Reference database</b>	Crystallographic Open Database

Table 1 Measurement parameters for LaB<sub>6</sub> sample

The diffractometer settings are constant for the **planet**. The operator can choose the measurement range and the integration time for optimal results and ease-of-use.

---

## Results

Figure 1 shows the complete diffraction pattern we recorded from the as received LaB<sub>6</sub> specimen.




Figure 1 Diffraction pattern from commercial LaB<sub>6</sub>

To assess the attainable resolution of the **planet** we performed careful profile fitting on each observed reflection from LaB<sub>6</sub> and we reported the FWHM. The results of this exercise are shown in Table 2.

Table 2 Full Width at Half Maximum for LaB<sub>6</sub> reflections

Scattering angle ( $^{\circ}2\theta$ )	FWHM ( $^{\circ}2\theta$ )
21.357	0.198
30.384	0.109
37.441	0.092
43.506	0.084
48.956	0.070
53.987	0.079
63.217	0.083
67.546	0.084
71.744	0.091
75.842	0.089
83.843	0.102
87.789	0.087
99.639	0.083

As can be seen the FWHM is below 0.11° for all reflections except the first one (0.20°), where the beam path is significantly shorter. The length of the beam and the focus dimensions determine the attainable

resolution for any reflection. This confirms that the attainable resolution for all reflections is best in class, for any portable XRD in the market.

#### Comparison with a laboratory-based system

Because we used commercial LaB<sub>6</sub>, we had to validate that this material is suitable for resolution measurements. To this end we measured both our commercial sample and a standard material (NIST 660a) on a laboratory-based system. These measurements are given by the Red (commercial) and Black (standard) curve in Figure 2.




Figure 2 FWHM measurements on laboratory system and planet

In the same graph we plotted the FWHM values we measured on the **planet** (Blue curve). We determined the FWHM by profile fitting in “Match!” from Crystal Impact. Please note that, although we measured the commercial sample with the **planet**, the FWHM values of all but the first reflections are even below the ones measured on the standard sample with a normal laboratory-based diffractometer! Moreover, we see the best resolution between 40 and 65°2θ, where it matters most, because generally the peak density is highest in that range.

## **Conclusion**

The attainable resolution of the planet compares well to the attainable resolution of a standard laboratory-based instrument. The attainable resolution of the planet is best in class among currently available portable X-ray diffractometers.

